

IOL - Injector - OVD Compatibility chart

Bi-Flex

		MEDJET MC ^{1.6}	MEDJET MB ^{1.8}	MEDJET MA ^{2.2}	MEDJET MX ²⁻⁴		Compatible OVDs
877FAB / 877FABY	DIOPTR	Not recommended		0.0D - +35.0D			Visco-MC (HPMC)
677AB / 677ABY	DIOPTR	Not recommended	-10.0D - +29.5D	-10.0D - +45.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677AD / 677ADY	DIOPTR	Not recommended	-10.0D - +29.5D	-10.0D - +45.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677TA / 677TAY	SPHERICAL EQUIVALENT	Not recommended		-10.0D - +35.0D	-10.0D - +30.0D	+30.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D	+1.0D - +10.0D	+1.00D - +4.5D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677M / 677MY	DIOPTR	Not recommended	+8.0D - +29.5D	+8.0D - +35.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677MT / 677MTY	SPHERICAL EQUIVALENT	Not recommended		+8.0D - +35.0D	+8.0D - +9.5D	+10.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D	+1.0D - +4.5D	+1.0D - +6.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)

Z-Flex

		MEDJET MC ^{1.6}	MEDJET MB ^{1.8}	MEDJET MA ^{2.2}	MEDJET MX ²⁻⁴		Compatible OVDs
860FAB / 860FABY	DIOPTR	Not recommended		0.0D - +35.0D			Visco-MC (HPMC)
690AB / 690ABY	DIOPTR	Not recommended	-10.0D - +29.5D	-10.0D - +45.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
690AD / 690ADY	DIOPTR	Not recommended	-10.0D - +29.5D	-10.0D - +45.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
690TA / 690TAY	SPHERICAL EQUIVALENT	Not recommended		-10.0D - +35.0D	-10.0D - +30.0D	+30.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D	+1.0D - +10.0D	+1.00D - +4.5D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
690M / 690MY	DIOPTR	Not recommended	+8.0D - +29.5D	+8.0D - +35.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
690MT / 690MTY	SPHERICAL EQUIVALENT	Not recommended		+8.0D - +35.0D	+8.0D - +9.5D	+10.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D	+1.0D - +4.5D	+1.0D - +6.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)

Q-Flex

		MEDJET MC ^{1.6}	MEDJET MB ^{1.8}	MEDJET MA ^{2.2}	MEDJET MX ²⁻⁴		Compatible OVDs
640AB / 640ABY	DIOPTR	Not recommended	0.0D - +29.5D	0.0D - +35.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
640AD / 640ADY	DIOPTR	Not recommended	0.0D - +29.5D	0.0D - +35.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
640MY	DIOPTR	Not recommended	+8.0D - +29.5D	+8.0D - +35.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)

Other

		MEDJET MC ^{1.6}	MEDJET MB ^{1.8}	MEDJET MA ^{2.2}	MEDJET MX ²⁻⁴		Compatible OVDs
611HPS	DIOPTR	Not recommended		-10.0D - +35.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
18ALY	DIOPTR	Not recommended	+10.0D - +30.0D				Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)

1stQ AddOn

		MEDJET MC ^{1.6}	MEDJET MB ^{1.8}	MEDJET MA ^{2.2}	MEDJET MX ²⁻⁴		Compatible OVDs	
A46R	DIOPTR	Not recommended		-10.0D - +10.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)	
A45RT	SPHERICAL EQUIVALENT	Not recommended		-10.0D - +10.0D	-10.0D - 3.25D	-3.0D - +8.0D	+8.25D - +10.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D	+1.0D - +10.0D	+1.0D - +11.0D	+1.0D - +10.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
A45RD2	DIOPTR	Not recommended		-7.0D - +6.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)	
A45DT	SPHERICAL EQUIVALENT	Not recommended		-3.0D - +3.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)	
	CYLINDER			+1.0D - +4.5D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)	
A45SML	DIOPTR	Not recommended		-4.0D - +4.0D			Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)	

Preloaded Injection Systems			
Bi-Flex			
POB-MA Preloaded Injection System			
877PA / 877PAY	DIOPTEP	-10.0D - +35.0D	Visco-MC (HPMC) BioVis1.6% (NaHa)
877PEY	DIOPTEP	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6% (NaHa)
PIL-MA Preloaded Injection System			
677P / 677PY	DIOPTEP	0.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677PMY	DIOPTEP	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
Medicel Accuject PRO 2.1P Preloaded Injection System			
677CMY	DIOPTEP	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)

Z-Flex			
POB-MA Preloaded Injection System			
860PA / 860PAY	DIOPTEP	-10.0D - +35.0D	Visco-MC (HPMC) BioVis1.6% (NaHa)
860PEY	DIOPTEP	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6% (NaHa)

Q-Flex			
PIL-MA Preloaded Injection System			
640P / 640PY	DIOPTEP	0.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
640PM	DIOPTEP	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)