

IOL - Injector - OVD Compatibility chart

Bi-Flex								
		MEDJET MC <sup>1,6</sup>	MEDJET MB <sup>1,8</sup>	MEDJET MA <sup>2,2</sup>	MEDJET MX <sup>2,4</sup>		Medicel Accuject 2.1-1P Top Loaded Injector	Compatible OVDs
877FAB / 877FABY	DIOPTER	Not recommended		0.0D - +35.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
877EBY	DIOPTER	Not recommended		+8.0D - +35.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677AB / 677ABY	DIOPTER	Not recommended	-10.0D - +29.5D	-10.0D - +45.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677AD / 677ADY	DIOPTER	Not recommended	-10.0D - +29.5D	-10.0D - +45.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677TA / 677TAY	SPHERICAL EQUIVALENT	Not recommended		-10.0D - +35.0D	-10.0D - +30.0D	+30.0D - +35.0D	-10.0D - 35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D	+1.0D - +10.0D	+1.0D - +4.5D		
677M / 677MY	DIOPTER	Not recommended	+8.0D - +29.5D	+8.0D - +35.0D			+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677MT / 677MTY	SPHERICAL EQUIVALENT	Not recommended		+8.0D - +35.0D	+8.0D - +9.5D	+10.0D - +35.0D	Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D	+1.0D - +4.5D	+1.0D - +6.0D		

Z-Flex								
		MEDJET MC <sup>1,6</sup>	MEDJET MB <sup>1,8</sup>	MEDJET MA <sup>2,2</sup>	MEDJET MX <sup>2,4</sup>		Medicel Accuject 2.1-1P Top Loaded Injector	Compatible OVDs
860FAB / 860FABY	DIOPTER	Not recommended		0.0D - +35.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
690AB / 690ABY	DIOPTER	Not recommended	-10.0D - +29.5D	-10.0D - +45.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
690AD / 690ADY	DIOPTER	Not recommended	-10.0D - +29.5D	-10.0D - +45.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
690TA / 690TAY	SPHERICAL EQUIVALENT	Not recommended		-10.0D - +35.0D	-10.0D - +30.0D	+30.0D - +35.0D	Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D	+1.0D - +10.0D	+1.0D - +4.5D		
690MY	DIOPTER	Not recommended	+8.0D - +29.5D	+8.0D - +35.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
690MTY	SPHERICAL EQUIVALENT	Not recommended		+8.0D - +35.0D	+8.0D - +9.5D	+10.0D - +35.0D	Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D	+1.0D - +4.5D	+1.0D - +6.0D		

Q-Flex								
		MEDJET MC <sup>1,6</sup>	MEDJET MB <sup>1,8</sup>	MEDJET MA <sup>2,2</sup>	MEDJET MX <sup>2,4</sup>		Medicel Accuject 2.1-1P Top Loaded Injector	Compatible OVDs
640AB / 640ABY	DIOPTER	Not recommended	0.0D - +29.5D	0.0D - +35.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
640AD / 640ADY	DIOPTER	Not recommended	0.0D - +29.5D	0.0D - +35.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
640MY	DIOPTER	Not recommended	+8.0D - +29.5D	+8.0D - +35.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)

Other								
		MEDJET MC <sup>1,6</sup>	MEDJET MB <sup>1,8</sup>	MEDJET MA <sup>2,2</sup>	MEDJET MX <sup>2,4</sup>		Medicel Accuject 2.1-1P Top Loaded Injector	Compatible OVDs
611HPS	DIOPTER	Not recommended		-10.0D - +35.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
18ALY	DIOPTER	Not recommended	+10.0D - +30.0D			Not recommended		Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)

1stQ AddOn								
		MEDJET MC <sup>1,6</sup>	MEDJET MB <sup>1,8</sup>	MEDJET MA <sup>2,2</sup>	MEDJET MX <sup>2,4</sup>		Medicel Accuject 2.1-1P Top Loaded Injector	Compatible OVDs
A46R	DIOPTER	Not recommended		-10.0D - +10.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
A45RT	SPHERICAL EQUIVALENT	Not recommended		-10.0D - +10.0D	-10.0D - 3.25D	-3.0D - +8.0D	Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D	+1.0D - +10.0D	+1.0D - +11.0D		
A45RD2	DIOPTER	Not recommended		-7.0D - +6.0D			-5.0D - +5.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
A45DT	SPHERICAL EQUIVALENT	Not recommended		-3.0D - +3.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
	CYLINDER			+1.0D - +4.5D				
A45SML	DIOPTER	Not recommended		-4.0D - +4.0D			Not recommended	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)

Preloaded Injection Systems			
<b>Bi-Flex</b>			
POB-MA Preloaded Injection System			
877PA / 877PAY	DIOPTER	-10.0D - +35.0D	Visco-MC (HPMC) BioVis1.6% (NaHa)
877PEY	DIOPTER	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6% (NaHa)
PIL-MA Preloaded Injection System			
677P / 677PY	DIOPTER	0.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677PMY	DIOPTER	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
Medicel Accuject PRO 2.1-1P Preloaded Injection System			
677CMY	DIOPTER	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677CMTY	DIOPTER	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
677CTA / 677CTAY	DIOPTER	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
<b>Z-Flex</b>			
POB-MA Preloaded Injection System			
860PA / 860PAY	DIOPTER	-10.0D - +35.0D	Visco-MC (HPMC) BioVis1.6% (NaHa)
860PEY	DIOPTER	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6% (NaHa)
<b>Q-Flex</b>			
Medicel Accuject PRO 1.8-1P Preloaded Injection System			
640CMY	DIOPTER	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
PIL-MA Preloaded Injection System			
640P / 640PY	DIOPTER	0.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)
640PM	DIOPTER	+8.0D - +35.0D	Visco-MC (HPMC) BioVis1.6%,1.8%,3.0% (NaHa)